

Hepatitis B (anti-HBc)

Overview

Clinical

Hepatitis B, caused by the Hepatitis B Virus (HBV), poses a significant worldwide health threat. This viral infection leads to various liver diseases, encompassing acute and chronic hepatitis, cirrhosis, and primary liver cancer. Transmission occurs through contact with infected bodily fluids, including blood, saliva, vaginal secretions, and semen. Additionally, maternal transmission to newborns is a recognized route of infection¹.

In clinical settings the measurement of Antibodies to Hepatitis B core antigen (anti-HBc) is used as an aid to assess both acute and chronic Hepatitis B infection². People who have immunity to hepatitis B from a vaccine do not develop anti-HBc³.

Epidemiology


It is estimated that 30.4 million people (10.5% of all people estimated to be living with hepatitis B) were aware of their infection, while 6.6 million (22%) of the people diagnosed were on treatment¹.

According to last US CDC recommendations, it is recommended to test all adults aged 18 years and older for HBsAg , anti-HBs and anti-HBc at least once in their lifetime⁴.

Anti-HBc CLIA

Assay Scheme

Qualitative measurement of total antibodies to Hepatitis B Core antigen (anti-HBc) in human serum or plasma



1 Hepatitis B. WHO Factsheet. July 2023. Accessed March 2024. <https://www.who.int/news-room/fact-sheets/detail/hepatitis-b>

2. Guidelines on Hepatitis B and C testing. WHO. February 2017

3. Interpretation of Hepatitis B Serologic Test Results. CDC . July 2023. Accessed March 2024.

<https://www.cdc.gov/hepatitis/hbv/interpretationOfHepBSerologicResults.htm>

4. Screening and Testing for Hepatitis B Virus Infection: CDC Recommendations — United States, 2023. Accessed March 2024.

https://www.cdc.gov/mmwr/volumes/72/rr/rr7201a1.htm?s_cid=rr7201a1_w

Evaluation of anti-HBc CLIA vs reference assay

| Anti-HBc CLIA Assay | Specimens Tested | Specificity |
|---------------------|------------------|-------------|
| Blood Donors | 5097 | 99.9% |
| NEG | 210 | 100% |

Table 2: Specificity assessment was based upon testing unselected blood donor serum samples, including first time donors and hospitalized patients. Specificity was calculated with the results shown in the table above

| Anti-HBc CLIA Assay | Specimens Tested | Specificity |
|---------------------|------------------|-------------|
| EIA method | 400 | 100% |
| CLIA method | 205 | 100% |

Table 3: A panel of samples from different sources with verified anti-HBc positivity was tested with our anti-HBc CLIA assay in comparison with commercially available EIA and CLIA anti-HBc methods. The sensitivity results are shown in the table above.

Cross-reactivity Test with anti-HBc CLIA

| Cross-reactivity | |
|--|-----------|
| Cross-reactant type | Agreement |
| anti-Toxoplasma (<i>Toxoplasma gondii</i>) | 4/4 |
| anti-Rubella | 10/10 |
| anti-HIV (Human Immunodeficiency Virus) | 7/7 |
| anti-HSV (Herpes Simplex Virus) | 5/5 |
| anti-HAV IgG (Hepatitis A Virus) | 6/6 |
| anti-HAV IgM (Hepatitis A Virus) | 5/5 |
| anti-HDV IgG (Hepatitis delta Virus) | 5/5 |
| anti-HEV (Hepatitis E Virus) | 4/4 |
| anti-EBV (Epstein-Barr Virus) | 4/4 |
| Anti-CMV (Cytomegalovirus) | 6/6 |
| anti-VZV (Varicella Zoster Virus) | 5/5 |

Table 4. Cross-reactant sample testing. 61 specimens with potential cross-reactivity with the anti-HBc CLIA assay were tested against commercially available anti-HBc assay. Table above is showing the agreement between methods

Werfen's Biomaterial offering

Recombinant HBc Antigen (ref 3000-5100 / 3000-5102)

Storage: -20°C

Source: *Escherichia coli*

Storage buffer: PBS, pH 7.4

Protein concentration: 0.2 – 1.0 mg/mL

Preservative: None

The content within this brochure is provided for informational purposes.

Contact oem@werfen.com for further technical information and product availability